

# Capnostream® 20p Bedside Capnography Monitor

The only integrated, complete picture of oxygenation and ventilation<sup>1</sup>



**COVIDIEN**

*positive results for life™*



## Smart Capnography™ and Pulse Oximetry Technology

Smart Capnography is a suite of algorithms proven to reduce alarms and simplify the use of capnography monitoring. It includes Smart Breath Detection™, Smart Alarm for Respiratory Analysis™ (SARA) and Nellcor™ SatSeconds to support alarm management. Additionally, Smart Capnography offers algorithms that provide workflow solutions, including the Integrated Pulmonary Index™ (IPI) and the new Apnea Sat-Alert algorithm. With all features combined, clinicians have access to a broader range of respiratory status information and alarm management technology to help enhance patient safety and improve clinical efficiency.<sup>2-3</sup>

The Capnostream® 20p patient monitor is built on a legacy of proven performance. For nearly two decades, clinicians have relied on Microstream®-enabled capnography monitoring for an accurate, continuous view of ventilation adequacy on intubated and non-intubated patients, from neonate to adult.

The 'p' in Capnostream® 20p stands for 'plus', and integrates an expanded electronics platform for enhanced Smart Capnography algorithms. Apnea Sat-Alert, the first algorithm to be supported on the new 'plus' platform, is designed to alert caregivers to repetitive patterns of apneas and oxygen desaturation.


The platform also supports the integration of the microMediCO<sub>2</sub>™ Module, the core component of Microstream etCO<sub>2</sub> measurement technology. Smaller and requiring less power to operate than previous modules, the microMediCO<sub>2</sub> Module can be integrated into an even wider range of OEM host monitor configurations.



# Patient Information + Analysis + Efficiency = Smart

## Apnea-Sat Alert (ASA)\*

ASA tracks and reports apneas per hour (A/hr), and the oxygen desaturation index (ODI) indicates the ‘dips’ in SpO<sub>2</sub> (number of times the SpO<sub>2</sub> value dropped 4 percent or more from baseline and returned to baseline in 240 seconds or less). A visual alert appears if the A/hr exceed a preset threshold over user selectable time periods of 2, 4, 8 or 12 hours. A/hr and ODI are displayed in real time on the monitor home screen, and the data is available in trend reports, in print outs and through data export.



ASA trend screen displays an analysis of A/hr by length of the apneas.

\*Apnea-Sat Alert is FDA cleared for adults 22 years of age or older.

## Smart Breath Detection™ (SBD)


Proprietary filter and pattern recognition algorithm screens out low-amplitude “non-breath” etCO<sub>2</sub> excursions like snoring, talking or crying, to provide a more accurate respiratory rate.

## Smart Alarm for Respiratory Analysis™ (SARA)

Reducing distractions from clinically insignificant alarms helps preserve caregiver alarm vigilance, leading to improved patient safety.<sup>4</sup> Functioning in combination with the SBD algorithm, the SARA algorithm is proven to reduce clinically insignificant respiratory rate alarms by dynamically adjusting the respiratory rate averaging algorithm during periods of breath-to-breath cycle variability.<sup>5</sup>

## Integrated Pulmonary Index™ (IPI)

IPI incorporates four real-time respiratory measurements into a single number, displayed on a scale from 1 to 10, representing an inclusive respiratory profile.



Helpful in busy clinical environments, IPI provides a simple and comprehensive indication of respiratory status and trends, promoting early awareness of changes to a patient’s breathing.

| IPI | Patient Status                                  |
|-----|---|
| 10  | Normal  |
| 8-9 | Within normal range                             |
| 7   | Close to normal range; requires attention       |
| 5-6 | Requires attention and may require intervention |
| 3-4 | Requires intervention                           |
| 1-2 | Requires immediate intervention                 |

## Nellcor™ SatSeconds

SatSeconds analyzes desaturation events by multiplying their duration, in seconds, by the number of percentage points the patient drops below the SpO<sub>2</sub> alarm limit. As a safety precaution, when three or more SpO<sub>2</sub> alarm violations occur within 60 seconds, an alarm will sound even if the SatSeconds limit has not been reached.

## LoSat Expanded Accuracy

Nellcor™ adhesive sensors with OxiMax™ technology and LoSat Expanded Accuracy range give clinicians the ability to assess patients with critically low saturation readings (60%-80%), such as infants with congenital heart disease. Accurate readings at low saturation levels expand options for clinicians to effectively monitor pulse oximetry.<sup>6</sup>



# Remote Alarm Management Connectivity Capabilities

The Capnostream® 20p patient monitor platform supports a range of remote alarm annunciation connectivity options to help clinicians monitor patients through nurse call systems, stand-alone central stations, wireless phones and pagers.

## Nurse Call System Alarm Annunciation

By attaching a cable to an auxiliary nurse call panel in the patient's room, the Capnostream 20p patient monitor alarm will be sent to the institution's nurse call system. In more advanced nurse call systems, alarms can be annotated on pagers and wireless phones.



## Vital Sync™ Virtual Patient Monitoring Platform 2.4 (VPMP)

Vital Sync VPMP by Covidien offers Electronic Medical Record (EMR) connectivity and remote continuous patient monitoring. Clinicians can remotely view patient information from ventilators, capnography monitors and pulse oximeters on any web-enabled devices. Patient information is sent to the EMR, Clinical Information System, and alarm forwarding systems.

## Nuvon VEGA™\*, iSirona™\* and Philips VueLink Systems

The Capnostream 20p patient monitor sends all measured patient information, including IPI data, to Nuvon VEGA™\* systems and iSirona™\* systems for electronic charting. Physiologic data automatically flows to the clinical information system, at the desired interval. When connected to Philips VueLink interfacing module†, parameter information and alarm status of the Capnostream 20p is automatically integrated. Alarms are annunciated on the Philips patient monitor and connected information center.

\*Visual only

# Summary of Capnostream® 20p Patient Monitor Remote Alarm Management and Connectivity Options

| FUNCTIONALITY  | OPTIONS                                    |                     |                              |                  |  |                      |
|--|--|---------------------|------------------------------|------------------|--|----------------------|
|  | Basic Nurse Call                           | Advanced Nurse Call | Nuvon VEGA™*                 | iSirona™*        | Philips VueLink                              | Vital Sync™ VPMP     |
|  | Capnostream 20p Patient Monitor Connection |                     |                              |                  |  |                      |
|  | Nurse Call Cable                           | Nurse Call Cable    | Wireless or Ethernet Adapter | Ethernet Adapter | Capnostream 20 Monitor/Philips VueLink Cable | Wireless or Ethernet |
| Basic alarm annunciation on nurse call system                      | ●  | ●                   |                              |                  |  |                      |
| Alarm annunciation on existing hospital wireless pagers and phones |  | ●                   |                              |                  |  |                      |
| IPI display  |  |                     |                              |                  | ●  | ●                    |
| Alarm differentiation by type of alarm                             |  |                     |                              |                  | ●  | ●                    |
| Central station view with parameters and alarms                    |  |                     |                              |                  | ●†   | ●                    |
| Wireless paging (stand-alone, departmental system)                 |  |                     |                              |                  |  | ●                    |
| Paging to existing hospital paging system                          |  |                     |                              |                  | ●  | ●                    |
| Vitals and ADT HL7 connectivity                                    |  |                     | ●                            | ●                | ●  | ●                    |

†Visual only

# Specifications

| POWER SUPPLY  |  |
|---|--|
| Input Voltage   | 100-240VAC, 50/60Hz  |
| Fuses   | Two F3.15A 250 Volt  |
| Input Power   | 90 VA  |
| BATTERY   |  |
| Battery Type  | 14.8V, 4Ah Lithium-Ion   |
| Battery Operation   | 2.5h (without thermal recorder)  |
| Battery Charging Time   | 100% in 12h  |
| CONTROLS  |  |
| Front Panel   | 1 Switch for monitor On/Off control<br>4 specific function keys<br>1 optical encoder with switch   |
| DISPLAY   |  |
| Screen  | 162mm (6.4in) Color TFT Display<br>Pixel Pitch: 0.204 (H) x 0.204(V) mm (0.008in)<br>Active Display Area: 130.56 (H) x 97.92 (V) mm (5.14in x 3.86in)<br>Resolution 640 x 480 pixels<br>Viewing angle (vertical) 110°<br>Viewing angle (horizontal) 140° |
| Trace Speed   | 3.0, 6.3, 12.5 and 25 mm/sec   |
| Waveform sampling rate  | 75.7 samples/sec for SpO <sub>2</sub> (fixed)<br>20 samples/sec for Capnography (fixed)  |
| Trend Storage   | 8640 point storage <ul style="list-style-type: none"> <li>• 2h at 5s resolution</li> <li>• 24h at 10s resolution</li> <li>• 72h at 30s resolution</li> </ul>   |
| Trend Display   | Graphical Display: 2h, 6h, 12h views<br>Tabular Display: 60 min, 15 min, 3 min, 1.5 min, and minimum resolution (minimum resolution settable to 5, 10, or 30 seconds)  |
| MICROSTREAM® CAPNOGRAPHY                                      |  |
| CO <sub>2</sub> Units   | mmHg or kPa or Vol%  |
| CO <sub>2</sub> , etCO <sub>2</sub> , FiCO <sub>2</sub> Range | 0-150 mmHg   |
| CO <sub>2</sub> Waveform Resolution                           | 0.1 mmHg   |
| EtCO <sub>2</sub> , FiCO <sub>2</sub> Resolution              | 1 mmHg   |
| CO <sub>2</sub> Accuracy                                      | 0-38 mmHg: ± 2 mmHg<br>39-150 mmHg: ± (5% of reading + 0.08% for every 1 mmHg above 38 mmHg)   |
| Respiration Rate Range  | 0-150 bpm  |
| Respiration Rate Accuracy                                     | 0-70 bpm: ±1 bpm<br>71-120 bpm: ±2 bpm<br>121-150 bpm: ±3 bpm  |

| CO <sub>2</sub> Alarms                      | No breath, EtCO <sub>2</sub> high, EtCO <sub>2</sub> low, RR high, RR low, IPI low (IPI also requires pulse oximetry information) |
|---|---|
| Flow Rate                                   | 50 (42.5 ≤ flow ≤ 65) ml/min, flow measured by volume   |
| Waveform Sampling                           | 20 samples/s  |
| Response Time                               | 2.95 s (typical)  |
| Initialization Time                         | 40 s (typical)  |
| Calibration Interval                        | Initially calibrate after 1,200 operating hours, then once a year or after 4,000 operating hours, whichever comes first           |
| NELLCOR™ OXIMAX™ PULSE OXIMETRY             |   |
| SpO <sub>2</sub> Measurement Range          | 1% to 100%  |
| SpO <sub>2</sub> Accuracy: Adult - Neonate  |   |
| Saturation (% SpO <sub>2</sub> ± 1 SD)      |   |
| 70% to 100% ± 2 digits; ± 3 digits (Motion) |   |
| 60% to 80% ± 3 digits                       |   |
| Low perfusion                               | 70% to 100% ± 2 digits  |
| Pulse Rate                                  | 20 to 250 bpm ± 3 digits  |
| Low perfusion                               | 20 to 250 bpm ± 3 digits  |
| Alarms                                      | Adjustable Alarm Limits<br>SpO <sub>2</sub> high, SpO <sub>2</sub> low, Pulse Rate high, Pulse Rate low                           |
| Sat Sec Range                               | 10, 25, 50, 100   |
| ALARMS                                      |   |
| High Priority Patient Warning Alarms        | Flashing Red LED<br>Flashing Red Numeric<br>High Priority Alarm beep pattern<br>Alarm Indication on Screen<br>Nurse Call          |
| Patient Caution Alarms                      | Flashing Yellow LED<br>Flashing Yellow Numeric  |
| Medium Priority Alarm                       | Flashing Yellow LED<br>Triple beep every thirty seconds<br>Alarm Indication on Screen<br>Nurse Call                               |
| Advisories                                  | Beep once<br>Advisory Indication on Screen  |
| Silent Advisories                           | Advisory Indication on Screen   |
| Alarm Volume Control                        | 5 steps   |
| Temporary Alarm Silence                     | All audible alarms silenced for 2 minutes   |

# Summary of Capnostream® 20p Patient Monitor Features and Enhancements

| NEW FEATURES   | INTENDED BENEFITS  |
|--|--|
| Apnea-Sat Alert  | Provides summary screens and USB/printer reports that indicate when apnea and oxygenation events occur. Offers apnea and oxygenation indexes for selected monitoring epochs.   |
| Alarm limit default settings   | The factory default alarm limit levels were adjusted in response to findings from a clinical user survey.  |
| Alarm limit electronic export and printout                                   | Urgent alarm limits can be exported to a USB Flash Memory Drive or printed out for documentation purposes.   |
| Instant demo mode  | Immediate activation facilitates sales demonstrations.   |
| "Parameter Standby" mode   | Suspends monitoring until a valid physiological signal is detected.<br>Designed to improve workflow by allowing the caregiver to suspend monitoring and alarms because a FilterLine® sampling line or SpO <sub>2</sub> sensor disconnection as been acknowledged and authorized. |
| "FilterLine® Disconnected" and "SpO <sub>2</sub> Sensor Disconnected" alarms | Using the "nurse call" output, the alarm alerts remote caregivers that a patient is no longer being monitored.   |

| ENHANCED FEATURE  | INTENDED BENEFITS  |
|---|--|
| "Trend Printout" updated to include full memory and resolution as displayed on "Tab Trend" screen | All trend data is printed in the display interval format, permitting complete documentation of procedural sedation cases.  |
| Permanent alarm silence on/off toggle   | If Permanent Alarm Silence Mode is enabled under institutional defaults, the user can now activate alarms by pressing the Alarm Silence Key. This mode may be desirable for highly supervised procedural sedation cases. |

1. Maddox RR, Oglesby H, Williams CK, Fields M, Danello S., Continuous respiratory monitoring and a "smart" infusion system improve safety of patient-controlled analgesia in the postoperative period. [http://www.ahrq.gov/downloads/pub/advances2/vol14/Advances-Maddox\\_111.pdf](http://www.ahrq.gov/downloads/pub/advances2/vol14/Advances-Maddox_111.pdf)
2. ECRI Institute. The Hazards of Alarm Overload: Keeping Excessive Physiologic Monitoring Alarms from Impeding Care. ECRI Guidance Article, March 2007.
3. Hockman S, Glembot T, Niebel K. Comparison of capnography derived respiratory rate alarm frequency using the SARA algorithm versus an established nonadaptive respiratory rate alarm management algorithm in bariatric surgical patients. *Resp Care* (Open Forum Abstracts). 2009;12.
4. The hazards of alarm overload. Keeping excessive physiologic monitoring alarms from impeding care. *Health Devices*. 2007;36(3):73-83.
5. Hockman S, Glembot T, Niebel K. Comparison of capnography derived respiratory rate alarm frequency using the SARA algorithm versus an established non-adaptive respiratory rate alarm management algorithm in bariatric surgical patients. [Open forum abstracts] *Resp Care*. 2009;12.
6. Bebout DE, Mannheimer PD, Wun C-C. Site-dependent differences in the time to detect changes in saturation during low perfusion. *Crit Care Med*. 2001;29(12):A115. [Abstract]



COVIDIEN, COVIDIEN with logo, Covidien logo and *positive results for life* are U.S. and internationally registered trademarks of Covidien AG. <sup>TM</sup>\* is a trademark of its respective owner. Other brands are trademarks of a Covidien company. ©2013 Covidien. 13-PM-0114(1)

6135 GUNBARREL AVENUE  
BOULDER, CO 80301  
800-635-5267

COVIDIEN.COM/RMS