

Palindrome™
Precision Chronic Dialysis
Catheter Family



MAHURKAR™
Acute Dialysis
Catheter Family



Argyle™
Peritoneal Dialysis
Catheter Family



**DIALYSIS
IS CRITICAL
CONFIDENCE
IS ESSENTIAL**

Dialysis Product Catalogue

Acute Vascular Access 04–08

MAHURKAR™ Straight Extension Dual Lumen Catheters	04
MAHURKAR™ Straight Extension Dual Lumen Catheter Kits	04
MAHURKAR™ Curved Extension Dual Lumen Catheters	05
MAHURKAR™ Straight Extension Dual Lumen Catheter Kits	05
MAHURKAR™ 13.5 Fr High Flow Catheter Kits	06
MAHURKAR™ 13.5 Fr High Flow Straight Extension Dual Lumen Catheter Kits	06
MAHURKAR™ 13.5 Fr High Flow Curved Extension Dual Lumen Catheter Kits	06
MAHURKAR™ 12 Fr High Pressure Curved Extension Triple Lumen Catheters	07
MAHURKAR™ 12 Fr High Pressure Curved Extension Triple Lumen Catheter Kits	07
Argyle™ 8 Fr Single Lumen Straight Extension Catheter	08
Argyle™ 8 Fr Single Lumen Catheter Kits	08

Chronic Vascular Access 09–16

Palindrome™ Symmetric Tip Dual Lumen Catheter Kits	09
Palindrome™ SI – Silver Ion Antimicrobial Dialysis Catheter Kits	10
Palindrome™ H – Heparin Coated Dialysis Catheter Kits	11
Palindrome™ HSI – Heparin Coated and Silver Ion Antimicrobial Dialysis Catheter Kits	12
Palindrome™ RT – Reverse Tunneled Dialysis Catheter Kits	13
Palindrome™ RT – Reverse Tunneled Dialysis Repair Catheter Kits	13
MAHURKAR™ 14.5 Fr Dual Lumen Catheter Kits	14
Palindrome™ and MAHURKAR™ Haemodialysis Repair Kits	14
Permcath™ Dual Lumen Catheters	15
Permcath™ Dual Lumen Catheter Kits	15
Permcath™ Dual Lumen Catheter Repair Kits	16
MAHURKAR™ 13.5 Fr Silicone Cuffed Catheter Kits	16

Peritoneal Access 17–21

Argyle™ Swan Neck Curl Cath Peritoneal Catheters.....	17
Argyle™ Swan Neck Curl Cath Missouri Peritoneal Catheters	17
Argyle™ Swan Neck Missouri Peritoneal Catheters.....	17
Argyle™ Moncrief-Popovich Swan Neck Curl Cath Peritoneal Catheter.....	17
Argyle™ Swan Neck Presternal Peritoneal Catheter	18
Argyle™ Swan Neck Tenckhoff Peritoneal Catheters.....	18
Argyle™ Swan Neck Oreopoulos-Zellerman Missouri Peritoneal Catheters.....	18
Oreopoulos-Zellerman Missouri Tenckhoff Peritoneal Catheter.....	18
Argyle™ Curl Cath Peritoneal Catheters – Two Cuffs.....	19
Argyle™ Curl Cath Peritoneal Catheters – One Cuff	19
Argyle™ Tenckhoff Peritoneal Catheters – Two Cuffs	19
Argyle™ Tenckhoff Peritoneal Catheters – One Cuff	19
Argyle™ Swan Neck Curl Cath Peritoneal Catheter Kits.....	20
Argyle™ Curl Cath and Tenckhoff Catheter Kits.....	20
Argyle™ Peri-Patch Repair Kit	21
Argyle™ Pediatric Peritoneal Catheters.....	21

Argyle™ Catheter Accessories 22–25

Argyle™ Pull-Apart Introducer Set with PTFE.....	22
Argyle™ Pull-Apart Introducer Sets.....	22
8 Fr/10 Fr Dilator.....	22
J/Straight Stainless Steel Guidewire	22
Argyle™ Y-Adapter	23
Argyle™ Introducer Needle	23
Argyle™ Injection Sealing Caps.....	23
Argyle™ Faller Tunneling Stylet/Bladed Trocar	23
Argyle™ Peritoneal Tunneling Stylet.....	24
Argyle™ Peritoneal Catheter Straightening Stylet.....	24
Argyle™ 2 Part Titanium Connector.....	24
Argyle™ Titanium Catheter Extender.....	24
Argyle™ Beta-Cap™ Adapter.....	24
Argyle™ Beta-Cap™ Clamp	25
Argyle™ Catheter Stencil – Swan Neck Tenckhoff.....	25
Argyle™ Catheter Stencil – Swan Neck Missouri.....	25
Medical Grade Silicone Adhesive Type A.....	25

Index 26

Acute Vascular Access



MAHURKAR™ Straight Extension Dual Lumen Catheters

- Flexible polyurethane material facilitates ease of insertion yet softens at body temperature
- Radiopaque polyurethane with a radiopaque tip for quick visualisation
- Patented Double-D™ lumen design optimises flow with lower arterial pressure
- Laser-cut side slots (11.5 Fr) designed to reduce the likelihood for positional occlusion and clot formation by minimising debris attachment
- Multiple configurations and sizes accommodate various insertion sites
- Soft, atraumatic green tip designed to reduce vessel trauma during insertion
- Clear silicone extensions provide kink resistance and view of blood
- Catheter extensions imprinted with name, length and priming volumes
- Durable Ultem™ adapters resist cracking
- Transparent, rotatable suture wing facilitates easy viewing of insertion site
- Sterile

SIDE HOLE CODES	SIDE SLOTS CODES	CATHETER SIZE (Fr)	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8817142005		10	12	0.8	0.9	5
8817145005		10	15	0.9	1.0	5
	8813816005	11.5	13.5	1.0	1.1	5
	8830414001	11.5	16	1.1	1.2	5
	8813794005	11.5	19.5	1.2	1.3	5
	8831662001	11.5	24	1.4	1.5	5



MAHURKAR™ Straight Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) Catheter (8, 10 or 11.5 Fr)
- (1) 18 ga Introducer Needle
- (1) J/Straight Guidewire - (0.032"/8 Fr, 0.035"/10 Fr or 0.038"/11.5 Fr)
- (1) 8 Fr or 10 Fr Dilator
- (1) 12 Fr Dilator (11.5 Fr Kits Only)
- (2) Telfa™ Island Dressings
- (2) Sealing Caps

SIDE HOLE CODES	SIDE SLOTS CODES	CATHETER SIZE (Fr)	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8832539001		8	9	0.8	0.9	5
8832539002		8	12	0.9	0.9	5
8832539003		8	15	1.0	1.0	5
8817146001		10	15	0.9	1.0	5
	8813817005	11.5	13.5	1.0	1.1	5
	8830415001	11.5	16	1.1	1.2	5
	8813793009	11.5	19.5	1.2	1.3	5
	8831661001	11.5	24	1.4	1.5	5



MAHURKAR™ Curved Extension Dual Lumen Catheters

- Flexible polyurethane material facilitates ease of insertion yet softens at body temperature
- Radiopaque polyurethane with a radiopaque tip for quick visualisation
- Patented curved extensions for jugular placement
- 180° curve bends away from patient's neck
- Patented Double-D™ lumen design optimises flow with lower arterial pressure
- Laser-cut side slots (11.5 Fr) designed to reduce the likelihood for positional occlusion and clot formation by minimising debris attachment
- Multiple configurations and sizes accommodate various insertion sites
- Soft, atraumatic green tip designed to reduce vessel trauma during insertion
- Clear silicone extensions provide kink resistance and view of blood
- Catheter extensions imprinted with name, length and priming volumes
- Durable Ultem™ adapters resist cracking
- Transparent, rotatable suture wing facilitates easy viewing of insertion site
- Sterile

CODE 11.5 Fr	CATHETER SIZE (Fr)	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8813816009	11.5	13.5	1.1	1.2	5
8830414002	11.5	16	1.2	1.3	5
8813794009	11.5	19.5	1.3	1.4	5



MAHURKAR™ Straight Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) 8 Fr, 10 Fr, or 11.5 Fr Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight Guidewire
- (1) 8 Fr Dilator (8 Fr catheter only)
- (1) 10 Fr Dilator (10 Fr and 11.5 Fr kits only)
- (1) 12 Fr Dilator (11.5 Fr catheter only)
- (2) Telfa™ Island Dressings
- (2) Sealing Caps

SIDE HOLE CODES	SIDE SLOTS CODES	CATHETER SIZE (Fr)	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8832539006		8	12	1.0	1.1	5
8832539007		8	15	1.1	1.2	5
8817143005		10	12	0.9	1.0	5
8817146007		10	15	1.0	1.1	5
8817149007		10	19.5	1.1	1.2	5
	8813817009	11.5	13.5	1.1	1.2	5
	8830415003	11.5	16	1.2	1.3	5
	8813793013	11.5	19.5	1.3	1.4	5

MAHURKAR™ 13.5 Fr High Flow Catheter Kits

- Flexible polyurethane material facilitates ease of insertion yet softens at body temperature
- Radiopaque polyurethane with a radiopaque tip for quick visualisation
- Larger Double-D™ lumen design optimises flow with lower arterial pressure
- Laser-cut side slots designed to reduce the likelihood for positional occlusion and clot formation by minimising debris attachment
- Soft, atraumatic green tip designed to reduce vessel trauma during insertion
- Clear silicone extensions provide kink resistance and view of blood
- Catheter extensions imprinted with name, length and priming volumes
- Durable Ultem™ adapters resist cracking
- Transparent, rotatable suture wing facilitates easy viewing of insertion site
- Sterile



MAHURKAR™ 13.5 Fr High Flow Straight Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) 13.5 Fr Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" × 70 cm Guidewire
- (1) 10 Fr Dilator
- (1) 14 Fr Dilator
- (2) Telfa™ Island Dressings
- (2) Sealing Caps

CODE 13.5 Fr	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8888135131	13.5	1.3	1.4	5
8888135161	16	1.4	1.5	5
8888135191	19.5	1.6	1.7	5
8888135241	24	1.8	1.9	5



MAHURKAR™ 13.5 Fr High Flow Curved Extension Dual Lumen Catheter Kits

Sterile Contents:

- (1) 13.5 Fr Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" × 70 cm Guidewire
- (1) 10 Fr Dilator
- (1) 14 Fr Dilator
- (2) Telfa™ Island Dressings
- (2) Injection Sealing Caps

CODE 13.5 Fr	IMPLANT LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8888135132	13.5	1.4	1.5	5
8888135162	16	1.5	1.6	5
8888135192	19.5	1.7	1.8	5
8888135242	24	1.9	2.0	5



Modified Double-D™ Design

MAHURKAR™ 12 Fr High Pressure Curved Extension Triple Lumen Catheters

- Indicated for dialysis, apheresis, infusion, high pressure injection & central venous pressure monitoring
- Flexible polyurethane material facilitates ease of insertion yet softens at body temperature
- Radiopaque polyurethane with a radiopaque tip for quick visualisation
- Laser-cut side slots designed to reduce the likelihood for positional occlusion and clot formation by minimising debris attachment
- Soft atraumatic tip designed to reduce vessel trauma during insertion
- 19 ga distinctive 3rd dedicated line
- Modified Double-D™ Lumen design ensures equal and consistent flow rates
- Imprinted priming volumes on catheter extensions quickly identify catheter locking volumes
- 180° curved extensions help stabilise the catheter and provide patient comfort
- Rotatable suture wings enable the tip to be rotated after implantation and stabilises the catheter
- Durable Ultem™ adapters resist cracking
- Sterile

LUMEN PRIMING VOLUMES (mL)

CODE 12 Fr	IMPLANT LENGTH (cm)	RED OR ARTERIAL	BLUE OR VENOUS	CLEAR OR MEDIAL	SHIP CASE
8888345504HP	13	1.4	1.4	0.4	5
8888345512HP	16	1.5	1.5	0.4	5
8888345520HP	20	1.6	1.7	0.5	5
8888345538HP	24	1.8	1.8	0.5	5

MAHURKAR™ 12 Fr High Pressure Curved Extension Triple Lumen Catheter Kits

Sterile Contents:

- (1) 12 Fr Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.035" × 70 cm Guidewire
- (1) 10 Fr Dilator
- (1) 12 Fr Dilator
- (2) Telfa™ Island Dressings
- (3) Sealing Caps

LUMEN PRIMING VOLUMES (mL)

CODE 12 Fr	IMPLANT LENGTH (cm)	RED OR ARTERIAL	BLUE OR VENOUS	CLEAR OR MEDIAL	SHIP CASE
Curved Extension					
8888345603HP	13	1.4	1.4	0.4	5
8888345611HP	16	1.5	1.5	0.4	5
8888345629HP	20	1.6	1.7	0.5	5
8888345637HP	24	1.8	1.8	0.5	5



Argyle™ 8 Fr Single Lumen Straight Extension Catheter

- Made of flexible radiopaque polyurethane
- Five outflow holes spiral near the tapered tip
- Rotatable transparent suture wings
- Luer-lock adapter
- Extension imprinted with priming volumes
- Sterile

CODE 8 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	PRIMING VOLUME (mL)	SHIP CASE
8817232019	19.5	29.0	1.2	5



Argyle™ 8 Fr Single Lumen Catheter Kits

Sterile Contents:

- (1) 8 Fr Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.035" Guidewire
- (1) Sealing Cap

CODE 8 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	PRIMING VOLUME (mL)	SHIP CASE
-----------	---------------------	---------------------	---------------------	-----------

Straight Extension

8831173010	15	24.5	1.0	5
8831173011	19.5	29.0	1.2	5

Curved Extension

8831173012	15	24.5	1.2	5
------------	----	------	-----	---

Chronic Vascular Access



Palindrome™ Symmetric Tip Dual Lumen Catheter Kits

- Available in the standard and precision tip configuration with radiopaque tungsten marker band, enhanced laser-cut side slots and optimised tip geometry.

Sterile Contents:

- (1) 14.5 Fr Symmetric Tip Catheter
- (2) Tal VenaTrac™ Over-the-Wire Insertion Stylets (additional option)
- (1) 16 Fr Valved Pull-Apart Sheath/Dilator
- (1) Bifurcated Tunneler
- (1) 12 Fr Tissue Dilator
- (1) 14 Fr Tissue Dilator
- (2) Sealing Caps
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (2) Telfa™ Island Dressing
- (4) 4" × 4" Cotton Gauze Sponges

CODE 14.5 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
--------------	---------------------	---------------------	------------------------------	----------------------------	-----------

Standard tip

8888145014	19	36	1.6	1.6	5
8888145015	23	40	1.9	1.9	5
8888145016	28	45	2.1	2.1	5
8888145017	33	50	2.3	2.3	5
8888145018	55	72	3.1	3.1	5

Precision tip

8888145014P	19	36	1.6	1.6	5
8888145015P	23	40	1.9	1.9	5
8888145016P	28	45	2.1	2.1	5
8888145017P	33	50	2.3	2.3	5
8888146044P	44	61	2.6	2.6	5
8888145018P	55	72	3.1	3.1	5

Standard tip with Tal VenaTrac™ Insertion Stylets

8888145039	19	36	1.6	1.6	5
8888145040	23	40	1.9	1.9	5
8888145041	28	45	2.1	2.1	5
8888145042	33	50	2.3	2.3	5

Precision tip with Tal VenaTrac™ Insertion Stylets

8888145039P	19	36	1.6	1.6	5
8888145040P	23	40	1.9	1.9	5
8888145041P	28	45	2.1	2.1	5
8888145042P	33	50	2.3	2.3	5



Palindrome™ SI – Silver Ion Antimicrobial Dialysis Catheter Kits

- Antimicrobial silver ions work to reduce microbial colonisation on the catheter surface
- Available in the standard and precision tip configuration with radiopaque tungsten marker band, enhanced laser-cut side slots and optimised tip geometry

Sterile Contents:

- (1) 14.5 Fr Symmetric Tip Catheter
- (2) Tal VenaTrac™ Over-the-Wire Insertion Stylets
- (1) 16 Fr Valved Pull-Apart Sheath/Dilator
- (1) Bifurcated Tunneler
- (1) 12 Fr Tissue Dilator
- (1) 14 Fr Tissue Dilator
- (2) Sealing Caps
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (2) Telfa™ Island Dressing
- (4) 4" × 4" Cotton Gauze Sponges

CODE 14.5 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
Standard tip					
8888145062C	19	36	1.6	1.6	5
8888145063C	23	40	1.9	1.9	5
8888145064C	28	45	2.1	2.1	5
8888145065C	33	50	2.3	2.3	5
8888145066C*	55	72	3.1	3.1	5
Precision tip					
8888145062CP	19	36	1.6	1.6	5
8888145063CP	23	40	1.9	1.9	5
8888145064CP	28	45	2.1	2.1	5
8888145065CP	33	50	2.3	2.3	5
8888146144CP	44	61	2.6	2.6	5
8888145066CP*	55	72	3.1	3.1	5

* Does not include Tal VenaTrac™ insertion stylets.



Palindrome™ H – Heparin Coated Dialysis Catheter Kits

- Aim to decrease the likelihood of clot formation and inhibits fibrin sheath propagation with its non-eluting heparin coating
- Available in the standard and precision tip configuration with radiopaque tungsten marker band, enhanced laser-cut side slots and optimised tip geometry

Sterile Contents:

- (1) 14.5 Fr Symmetric Tip Catheter
- (2) Tal VenaTrac™ Over-the-Wire Insertion Stylets
- (1) 16 Fr Valved Pull-Apart Sheath/ Dilator
- (1) Bifurcated Tunneler
- (1) 12 Fr Tissue Dilator
- (1) 14 Fr Tissue Dilator
- (2) Sealing Caps
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (2) Telfa™ Island Dressing
- (4) 4" × 4" Cotton Gauze Sponges

CODE 14.5 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
Standard tip					
8888145043C	19	36	1.6	1.6	5
8888145044C	23	40	1.9	1.9	5
8888145045C	28	45	2.1	2.1	5
8888145046C	33	50	2.3	2.3	5
Precision tip					
8888145043CP	19	36	1.6	1.6	5
8888145044CP	23	40	1.9	1.9	5
8888145045CP	28	45	2.1	2.1	5
8888145046CP	33	50	2.3	2.3	5



Palindrome™ HSI – Heparin Coated and Silver Ion Antimicrobial Dialysis Catheter Kits

The Palindrome™ HSI – Heparin Coated and Silver Ion Antimicrobial Dialysis Catheter is the premier dialysis catheter, incorporating the non-eluting heparin coating to decrease the likelihood of clot formation and fibrin sheath dissemination along with the silver ion sleeve to reduce microbial colonisation on the catheter surface.

- Available in the standard and precision tip configuration with radiopaque tungsten marker band, enhanced laser-cut side slots and optimised tip geometry

Sterile Contents:

- (1) 14.5 Fr Symmetric Tip Catheter
- (2) Tal VenaTrac™ Over-the-Wire Insertion Stylets
- (1) 16 Fr Valved Pull-Apart Sheath/ Dilator
- (1) Bifurcated Tunneler
- (1) 12 Fr Tissue Dilator
- (1) 14 Fr Tissue Dilator
- (2) Sealing Caps
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (2) Telfa™ Island Dressing
- (4) 4" × 4" Cotton Gauze Sponges

CODE 14.5 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
Standard tip					
8888145057C	19	36	1.6	1.6	5
8888145048C	23	40	1.9	1.9	5
8888145049C	28	45	2.1	2.1	5
8888145050C	33	50	2.3	2.3	5
Precision tip					
8888145057CP	19	36	1.6	1.6	5
8888145048CP	23	40	1.9	1.9	5
8888145049CP	28	45	2.1	2.1	5
8888145050CP	33	50	2.3	2.3	5



Palindrome™ RT – Reverse Tunneled Dialysis Catheter Kits

- Allows for precise tip placement and ideal tunnel trajectory using the retrograde tunnel technique
- Available in the standard and precision tip configuration with radiopaque tungsten marker band, enhanced laser-cut side slots and optimised tip geometry

Sterile Contents:

- (1) Catheter/Valve Adapter Assembly
- (1) Hub/Back-end Assembly with Colored End Cap
- (1) Hub Snap Connector
- (1) Tunneler
- (1) Additional Tunneler Cap
- (1) Syringe
- (4) 4" x 4" Cotton Gauze Sponges
- (1) 12 Fr (4.0 mm) Dilator
- (1) 14 Fr (4.7 mm) Dilator
- (1) 16 Fr (5.3 mm) Valved Pull-Apart Sheath/Dilator
- (1) Telfa™ Island Dressing
- (2) Smooth Jawed Forceps
- (2) Sealing Caps
- (1) External Measuring Kit (includes 18 ga [1.2 mm] Introducer Needle, #11 scalpel, 0.038" [0.965 mm] J/straight guidewire)

CODE	INSERTION LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
------	-----------------------	---------------------	------------------------------	----------------------------	-----------

Standard tip

8888541019	19	39	1.5	1.5	5
8888541023	23	43	1.7	1.7	5
8888541028	28	48	1.9	1.9	5
8888541033	33	53	2.1	2.1	5
8888541044	44	64	2.5	2.5	5
8888541055	55	75	3.0	3.0	5

Precision tip

8888541019P	19	39	1.5	1.5	5
8888541023P	23	43	1.7	1.7	5
8888541028P	28	48	1.9	1.9	5
8888541033P	33	53	2.1	2.1	5
8888541044P	44	64	2.5	2.5	5
8888541055P	55	75	3.0	3.0	5



Palindrome™ RT – Reverse Tunneled Dialysis Repair Catheter Kits

Sterile Contents:

- (1) Hub/Back-end Assembly
- (1) Hub Snap Connector
- (1) Smooth Jawed Forceps
- (1) Sealing Caps
- (1) Drape
- (1) Ruler

CODE	IMPLANT LENGTH (cm)	SHIP CASE
8888541119	19	1
8888541123	23	1
8888541128	28	1
8888541133	33	1
8888541144	44	1
8888541155	55	1

Note: Ensure the catheter repair kit corresponds to the same implant length catheter.



MAHURKAR™ 14.5 Fr Dual Lumen Catheter Kits

Sterile Contents:

- (1) 14.5 Fr Catheter
- (1) 16 Fr Valved Pull-Apart Sheath/Dilator
- (1) Tunneling Stylet, 21 cm
- (1) 12 Fr Tissue Dilator
- (1) 14 Fr Tissue Dilator
- (2) Sealing Caps
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) #11 Scalpel
- (2) Telfa™ Island Dressing
- (4) 4" × 4" Cotton Gauze Sponges

CODE 14.5 Fr	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8888145251	19	39	1.6	1.7	5
8888145252	23	40	1.8	1.9	5
8888145253	28	45	2.0	2.1	5
8888145254	33	50	2.2	2.3	5



Palindrome™ and MAHURKAR™ Haemodialysis Repair Kits

Sterile Contents:

- (1) Arterial Repair Assembly
- (1) Venous Repair Assembly
- (2) Sealing Caps
- (2) Temporary Slide Clamps
- (1) Drape
- (1) Disposable Scissors
- (1) Measurement Guide
- (1) Priming Volume Label Sheet

CODE	DESCRIPTION	SHIP CASE
8888200001	Haemodialysis Repair Kit	1

Note: Repairs Palindrome™ chronic dialysis catheters (except Palindrome™ RT) and MAHURKAR™ acute and chronic catheters.



Permcath™ Dual Lumen Catheters

- Made of soft silicone rubber
- 2.5 cm separation between arterial lumen and venous tip
- May be implanted in the jugular, subclavian or femoral vein
- Flexible silicone suture wings
- Felt cuff anchors catheter subcutaneously
- Clear extensions imprinted with priming volumes
- Two Sealing Caps
- Pediatric option available

CODE	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8815132001	19	36	1.3	1.4	1
8817606001	23	40	1.4	1.5	1

Pediatric

8815543001	13	28	0.8	0.85	1
------------	----	----	-----	------	---



Permcath™ Dual Lumen Catheter Kits

Sterile Contents:

- (1) Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) 10 Fr Dilator (13 cm only)
- (1) 12 Fr Dilator (except 13 cm)
- (1) #11 Scalpel
- (1) Tunneling Stylet
- (1) Oval Pull-Apart Sheath/Dilator
- (4) 4" x 4" Cotton Gauze Sponges
- (2) Telfa™ Island Dressings
- (2) Sealing Caps

CODE	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8817748001	19	36	1.3	1.4	5
8817749001	23	40	1.4	1.5	5
8831692001	28	45	1.6	1.7	5

Pediatric

8834369001	13	28	0.8	0.85	5
------------	----	----	-----	------	---



Permcath™ Dual Lumen Catheter Repair Kits

- Replacement Section with Dual Extensions
- Luer-lock Adapters
- In-line Clamps
- Silicone Rubber Repair Sleeve

Sterile Contents:

- (1) Replacement Section with Dual Extensions
- (1) 1 mL Syringe with Blunted Needle
- (2) Sealing Caps

CODE	FOR PERMCATH CATHETER	SHIP CASE
8888680009	Permcath, 13 cm	1
8888169001	Permcath, 19 cm	1
8888169003	Permcath, 23 cm	1

Note: Repair Kits are not universal. They must be ordered for the specific size catheter. To be repaired, an indwelling catheter must have at least 3 cm external tubing at the exit site. Silicone Adhesive is packaged individually and must be ordered separately to ensure maximum shelf-life. See Page #25 for ordering information (Code #8810807001).



MAHURKAR™ 13.5 Fr Silicone Cuffed Catheter Kits

- Made of soft silicone rubber
- Radiopaque
- 2.5 cm separation between arterial lumen and venous tip
- May be implanted in the jugular, subclavian or femoral vein
- Flexible silicone suture wings
- Felt cuff anchors catheter subcutaneously
- Clear extensions imprinted with priming volumes

Sterile Contents:

- (1) Catheter
- (1) 18 ga Introducer Needle
- (1) J/Straight 0.038" Guidewire
- (1) 12 mL Syringe
- (1) 12 Fr Dilator
- (1) #11 Scalpel
- (1) Tunneling Stylet
- (1) Pull-Apart Sheath/Dilator
- (4) 4" x 4" Cotton Gauze Sponges
- (2) Telfa™ Island Dressings
- (2) Sealing Caps

CODE	IMPLANT LENGTH (cm)	OVERALL LENGTH (cm)	ARTERIAL PRIMING VOLUME (mL)	VENOUS PRIMING VOLUME (mL)	SHIP CASE
8833074001	19	36	1.4	1.5	5
8833074002	23	40	1.5	1.6	5
8833074003	28	45	1.7	1.8	5
8833074004	33	50	1.9	2.0	5

Peritoneal Access



Argyle™ Swan Neck Curl Cath Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888413807	Swan Neck Curl Cath Catheter 2 Cuff Left, 62.5 cm	1
8888413815	Swan Neck Curl Cath Catheter 2 Cuff Right, 62.5 cm	1



Argyle™ Swan Neck Curl Cath Missouri Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888413401	Swan Neck Curl Cath Missouri Catheter 2 Cuff Left, 62 cm	1
8888413419	Swan Neck Curl Cath Missouri Catheter 2 Cuff Right, 62 cm	1



Argyle™ Swan Neck Missouri Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888412601	Swan Neck Missouri Catheter 2 Cuff Left, 44.5 cm	1
8888412619	Swan Neck Missouri Catheter 2 Cuff Right, 44.5 cm	1



Argyle™ Moncrief-Popovich Swan Neck Curl Cath Peritoneal Catheter

CODE	CATHETER DESCRIPTION	SHIP CASE
8888414029	Swan Neck Moncrief-Popovich Swan Neck Curl Cath Catheter 2 Cuff, 62.5 cm	5



Argyle™ Swan Neck Presternal Peritoneal Catheter

CODE	CATHETER DESCRIPTION	SHIP CASE
8888414011	Swan Neck Presternal Catheter 2 Cuff, 112.8 cm (Swan Neck: 60.3 cm + 52.5 cm Curl)	1



Argyle™ Swan Neck Tenckhoff Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888412007	Swan Neck Tenckhoff Catheter 2 Cuff Left, 43 cm	1
8888412015	Swan Neck Tenckhoff Catheter 2 Cuff Right, 43 cm	1
8888412023	Swan Neck Tenckhoff Catheter 2 Cuff Left, 55 cm	1



Argyle™ Swan Neck Oreopoulos-Zellerman Missouri Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888413005	Swan Neck Oreopoulos-Zellerman Catheter, 2 Cuff Left, 44.5 cm	1
8888413013	Swan Neck Oreopoulos-Zellerman Catheter, 2 Cuff Right, 44.5 cm	1



Oreopoulos-Zellerman Missouri Tenckhoff Peritoneal Catheter

CODE	CATHETER DESCRIPTION	SHIP CASE
8888411009	Oreopoulos-Zellerman Missouri Tenckhoff Catheter, 1 Bead, 2 Disk, 2 Cuff, 41 cm	1



Argyle™ Curl Cath Peritoneal Catheters – Two Cuffs

CODE	CATHETER DESCRIPTION	SHIP CASE
8811313009	Curl Cath, Curled, 2 Mobile Cuff, 62 cm	1
8811313015	Curl Cath, Curled, 2 Cuff, 57 cm	1
8811313010	Curl Cath, Curled, 2 Cuff, 62 cm	1



Argyle™ Curl Cath Peritoneal Catheters – One Cuff

CODE	CATHETER DESCRIPTION	SHIP CASE
8811313014	Curl Cath, Curled, 57 cm	1
8811313013	Curl Cath, Curled, 60 cm	1
8888411710	Curl Cath, Curled, 62 cm	1



Argyle™ Tenckhoff Peritoneal Catheters – Two Cuffs

CODE	CATHETER DESCRIPTION	SHIP CASE
8810887003	Tenckhoff, 2 Mobile Cuff, 42 cm	1
8810888003	Tenckhoff, 2 Cuff, 42 cm	1
8810888012	Tenckhoff, Universal, 2 Cuff, 47 cm	1
8888411405	Tenckhoff, 2 Cuff, 41 cm	1
8888411447	Tenckhoff, 2 Cuff, 47 cm	1



Argyle™ Tenckhoff Peritoneal Catheters – One Cuff

CODE	CATHETER DESCRIPTION	SHIP CASE
8810889003	Tenckhoff, 1 Subcutaneous Cuff, 42 cm	1
8810889011	Tenckhoff, 1 Subcutaneous Cuff, 46 cm	1
8814843001	Tenckhoff, 1 Preperitoneal Cuff, 41 cm	1
8814843002	Tenckhoff, Universal 1 Preperitoneal Cuff, 46 cm	1
8888411439	Tenckhoff, 1 Cuff, 42 cm	1
8888411421	Tenckhoff, 1 Cuff, 47 cm	1



Argyle™ Swan Neck Curl Cath Peritoneal Catheter Kits

Sterile Contents:

- (1) Peritoneal Catheter
- (1) Beta-Cap Adapter, Cap and Clamp
- (1) 16 Fr Pull-Apart Sheath/Dilator
- (1) 18 ga Introducer Needle
- (1) J/Straight Guidewire
- (1) #11 Scalpel
- (1) 12 mL Syringe
- (1) Tunneling Stylet
- (6) 4" × 4" Gauze Sponges

CODE	CATHETER DESCRIPTION	SHIP CASE
8888413823	Swan Neck Curl Cath Kit, 2 Cuff, Left, 62.5 cm	1
8888413831	Swan Neck Curl Cath Kit, 2 Cuff, Right, 62.5 cm	1



Argyle™ Curl Cath Catheter Kit

Argyle™ Curl Cath and Tenckhoff Catheter Kits

Sterile Contents:

- (1) Peritoneal Catheter
- (1) Beta-Cap Adapter, Cap and Clamp
- (1) 16 Fr Pull-Apart Sheath/Dilator
- (1) 18 ga Introducer Needle
- (1) J/Straight Guidewire
- (1) #11 Scalpel
- (1) 12 mL Syringe
- (1) Tunneling Stylet
- (6) 4" × 4" Gauze Sponges

CODE	CATHETER DESCRIPTION	SHIP CASE
8817278007	Curl Cath, 2 Cuff, 57 cm	5
8817278006	Curl Cath, 2 Cuff, 62 cm	5
8817278010	Curl Cath, 1 Preperitoneal Cuff, 60 cm	5
8817278008	Tenckhoff, 2 Cuff, 42 cm	5
8817278001	Tenckhoff Universal, 1 Preperitoneal Cuff, 46 cm	5



Argyle™ Tenckhoff Catheter Kit



Argyle™ Peri-Patch Repair Kit

Sterile Contents:

- (1) 15 cm Silicone Rubber Catheter Extension (2.6 mm I.D., 4.9 mm O.D.) with Double-Barbed Connector
- (1) Beta-Cap Adapter, Cap and Clamp
- (1) Glue Mold

CODE	CATHETER DESCRIPTION	SHIP CASE
8888285001	Peri-Patch Repair Kit*	1

Note: Peri-Patch Repair Kits are universal for all Medtronic Peritoneal Dialysis Catheters. To be repaired, an indwelling catheter must have at least 3 cm external tubing at the exit site. Silicone Adhesive is packaged individually and must be ordered separately to ensure maximum shelf-life.

See Page #25 for ordering information (Code #8810807001).



Argyle™ Pediatric Peritoneal Catheters

CODE	CATHETER DESCRIPTION	SHIP CASE
8888413100	Pediatric Swan Neck Curl Cath Catheter, 2 Cuff, 38.9 cm	1
8888413101	Pediatric Swan Neck Curl Cath Catheter, 2 Cuff, 43 cm	1
8888413102	Adolescent Swan Neck Curl Cath Catheter, 2 Cuff, 59 cm	1
8888414813	Pediatric Swan Neck Curl Cath Catheter, 2 Cuff, Left, 42 cm	1
8888412643	Pediatric Swan Neck Missouri Catheter, 2 Cuff, Left, 38 cm	1
8888415307	Pediatric Swan Neck Curl Cath Missouri Catheter, 2 Cuff, Left, 42.5cm	1
8888410506	Pediatric Swan Neck Tenckhoff Catheter, 2 Cuff, Left, 37.5 cm	1
8888410514	Pediatric Swan Neck Tenckhoff Catheter, 2 Cuff, Right, 37.5 cm	1
8815677001	Pediatric Curl Cath Catheter, 1 Cuff, 39 cm	1
8888414201	Pediatric Tenckhoff, 2 Cuff, 30 cm	1
8888414219	Pediatric Tenckhoff, 1 Cuff, 30 cm	1
8812329002	Pediatric Tenckhoff Catheter, 2 Cuff, 31 cm	1
8810890014	Pediatric Tenckhoff Catheter, 2 Cuff, 32 cm	1
8810890003	Pediatric Tenckhoff Catheter, 2 Cuff, 37 cm	1
8812329001	Pediatric Tenckhoff Catheter, 1 Subcutaneous Cuff, 31 cm	1

Argyle™ Catheter Accessories



Argyle™ Pull-Apart Introducer Set with PTFE

Sterile Contents:

- (1) 16 Fr Pull-Apart Sheath/Dilator†
- (1) 18 ga Introducer Needle
- (1) 12 mL Syringe
- (1) J/Straight 0.038" Guidewire

CODE	DESCRIPTION	SHIP CASE
8815544010	16 Fr Pull-Apart Introducer Set†	5

† Made with PolyTetra FluoroEthylene (PTFE).




Argyle™ Pull-Apart Introducer Sets

Sterile Contents:

- (1) Pull-Apart Sheath/Dilator
- (1) 18 ga Introducer Needle
- (1) 12 mL Syringe
- (1) J/Straight 0.038" Guidewire
- (1) 12 Fr Dilator (Oval Only)

CODE	DESCRIPTION	SHIP CASE
8815544003	10 Fr Pull-Apart Introducer Set	5
8815544002	12 Fr Pull-Apart Introducer Set	5
8815544001	14 Fr Pull-Apart Introducer Set	5
8815544013	Oval Pull-Apart Introducer Set	5



8 Fr/10 Fr Dilator

- Assists catheter insertion by widening the vessel or surrounding tissue
- Luer lock hub
- Sterile

CODE	DESCRIPTION	SHIP CASE
8815072001	8 Fr Dilator, Length 15 cm	5
8815642001	10 Fr Dilator, Length 15 cm	5



J/Straight Stainless Steel Guidewire

- Sterile
- Guides dilator or catheter into a vein or peritoneal cavity
- Has 2 soft ends, 1 in a "J" configuration
- Tip straightener attached

CODE	DESCRIPTION	SHIP CASE
8813796001	J/Straight 0.038" × 70 cm Guidewire	10
8817231001	J/Straight 0.035" × 70 cm Guidewire	10



Argyle™ Y-Adapter

- Disposable
- Converts single extension to dual to allow for infusion of two compatible medications/solutions
- Sterile

CODE	DESCRIPTION	SHIP CASE
8818782001	Disposable Y-adapter	25



Argyle™ Introducer Needle

- Stainless steel
- Channels guidewire into a vein or peritoneal cavity
- Accommodates guidewires up to 0.038" diameter
- Luer lock hub
- Sterile

CODE	DESCRIPTION	SHIP CASE
8813795001	18 ga × 7 cm Introducer Needle	10



Argyle™ Injection Sealing Caps

- Luer lock fitting with a natural rubber stopper
- Provides a self-sealing injection site for administering medication or drawing blood
- Universal adapter
- May be used with both hemodialysis and peritoneal dialysis catheters
- Sterile

CODE	DESCRIPTION	SHIP CASE
8813791001	Injection Sealing Caps	50



Argyle™ Faller Tunneling Stylet/Bladed Trocar

- Sterile

CODE	DESCRIPTION	SHIP CASE
8888415679	Faller Tunneling Stylet/Bladed Trocar	1



- Argyle™ Peritoneal Tunneling Stylet
- Made of pliable stainless steel
- Can be bent to create customized, curved tunnel
- Sterile

CODE	DESCRIPTION	SHIP CASE
8817791001	5 Fr × 23 cm (9") Stylet	5



Argyle™ Peritoneal Catheter Straightening Stylet

- Stainless steel stylet
- Used to keep the catheter straight and maneuverable during insertion into the peritoneum
- Autoclavable
- Non-sterile

CODE	DESCRIPTION	SHIP CASE
8888415646	Insertion/Straightening Stylet, 44.5 cm	1
8888415661	Insertion/Straightening Stylet, 61 cm	1



Argyle™ 2 Part Titanium Connector

- Sterile

CODE	DESCRIPTION	SHIP CASE
8888415604	2 Part Titanium Connector	1



Argyle™ Titanium Catheter Extender

- Sterile

CODE	DESCRIPTION	SHIP CASE
8888415612	Adult Titanium Catheter Extender	1



Argyle™ Beta-Cap™ Adapter

- Made of Ultem™ resin
- For use with all Argyle™ Peritoneal Catheters
- Sterile

CODE	DESCRIPTION	SHIP CASE
8814661001	Beta-Cap Adapter	25



Argyle™ Beta-Cap™ Clamp

- For use with all Argyle™ Peritoneal Catheters
- Sterile

CODE	DESCRIPTION	SHIP CASE
8810805001	Beta-Cap Clamp	25



Argyle™ Catheter Stencil – Swan Neck Tenckhoff

- Non-sterile

CODE	DESCRIPTION	SHIP CASE
8888412304	Swan Neck Tenckhoff stencil, 2 Cuff, Left/Right	1



Argyle™ Catheter Stencil – Swan Neck Missouri

- Non-sterile

CODE	DESCRIPTION	SHIP CASE
8888412312	Swan Neck Missouri stencil, 2 Cuff, Left/Right	1



Medical Grade Silicone Adhesive Type A

- Minimum 72-hour curing required for catheter repair
- For use with all Argyle™ repair kits to extend or repair catheters
- Sterile

CODE	DESCRIPTION	SHIP CASE
8810807001	Silicone Adhesive, Type A	1

Index

CODE	PAGE	CODE	PAGE	CODE	PAGE
8810805001	25	8815544003	22	8832539003	04
8810807001	25	8815544010	22	8832539006	05
8810887003	19	8815544013	22	8832539007	05
8810888003	19	8815642001	22	8833074001	16
8810888012	19	8815677001	21	8833074002	16
8810889011	19	8817142005	04	8833074003	16
8810890003	21	8817143005	05	8833074004	16
8810890014	21	8817145005	04	8834369001	15
8811313009	19	8817146001	04	8888135131	06
8811313010	19	8817146007	05	8888135132	06
8811313013	19	8817149007	05	8888135161	06
8811313014	19	8817231001	22	8888135162	06
8811313015	19	8817232019	08	8888135191	06
8812329001	21	8817278001	20	8888135192	06
8812329002	21	8817278006	20	8888135241	06
8813791001	23	8817278007	20	8888135242	06
8813793009	04	8817278008	20	8888145014	09
8813793013	05	8817278010	20	8888145014P	09
8813794005	04	8817606001	15	8888145015	09
8813794009	05	8817748001	15	8888145015P	09
8813795001	23	8817749001	15	8888145016	09
8813796001	22	8817791001	24	8888145016P	09
8813816005	04	8818782001	23	8888145017	09
8813816009	05	8830414001	04	8888145017P	09
8813817005	04	8830414002	05	8888145018	09
8813817009	05	8830415001	04	8888145018P	09
8814661001	24	8830415003	05	8888145039	09
8814843001	19	8831173010	08	8888145039P	09
8814843002	19	8831173011	08	8888145040	09
8814843003	19	8831173012	08	8888145040P	09
8815072001	22	8831661001	04	8888145041	09
8815132001	15	8831662001	04	8888145041P	09
8815543001	15	8831692001	15	8888145042	09
8815544001	22	8832539001	04	8888145042P	09
8815544002	22	8832539002	04	8888145043C	11

Index

CODE	PAGE	CODE	PAGE	CODE	PAGE
8888145043CP	11	8888345504HP	07	8888414011	18
8888145044C	11	8888345512HP	07	8888414029	17
8888145044CP	11	8888345520HP	07	8888414201	21
8888145045C	11	8888345538HP	07	8888414219	21
8888145045CP	11	8888345603HP	07	8888414813	21
8888145046C	11	8888345611HP	07	8888415307	21
8888145046CP	11	8888345629HP	07	8888415604	24
8888145048C	12	8888345637HP	07	8888415612	24
8888145048CP	12	8888410506	21	8888415646	24
8888145049C	12	8888410514	21	8888415661	24
8888145049CP	12	8888411009	18	8888415679	23
8888145050C	12	8888411405	19	8888541019	13
8888145050CP	12	8888411421	19	8888541019P	13
8888145057C	12	8888411439	19	8888541023	13
8888145057CP	12	8888411447	19	8888541023P	13
8888145062C	10	8888411710	19	8888541028	13
8888145062CP	10	8888412007	18	8888541028P	13
8888145063C	10	8888412015	18	8888541033	13
8888145063CP	10	8888412023	18	8888541033P	13
8888145064C	10	8888412304	25	8888541044	13
8888145064CP	10	8888412312	25	8888541044P	13
8888145065C	10	8888412601	17	8888541055	13
8888145065CP	10	8888412619	17	8888541055P	13
8888145066C	10	8888412643	21	8888541119	13
8888145066CP	10	8888413005	18	8888541123	13
8888145251	14	8888413013	18	8888541128	13
8888145252	14	8888413100	21	8888541133	13
8888145253	14	8888413101	21	8888541144	13
8888145254	14	8888413102	21	8888541155	13
8888146044P	09	8888413401	17	8888680009	16
8888146144CP	10	8888413419	17		
8888169001	16	8888413807	17		
8888169003	16	8888413815	17		
8888200001	14	8888413823	20		
8888285001	21	8888413831	20		

Medtronic



Use scan app to read

For more information, please visit
medtronic.eu/product-catalog

IMPORTANT: Please refer to the package insert for complete instructions, contraindications, warnings and precautions.

© 2016 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic.

™* MAHURKAR, Double-D design and Double-D are trademarks of Sakharam D. Mahurkar, M.D. Used under license.

™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

16-eu-dialysis-product-catalogue-936586