


### Shiley™ PED, PDC, PLC and NEO

Cuffed and cuffless paediatric and neonatal tracheostomy tubes

#### PACK INCLUDES:

- All Shiley™ Neonatal and Paediatric Tracheostomy Tubes have a recessed area behind the connector to provide greater comfort.
- The Neonatal Tracheostomy Tube (NEO) features a lower flange angle to improve fit.
- The Long Paediatric Tracheostomy Tube (PDL) offers a longer cannula to fit larger paediatric patients.
- The economical Cuffed Paediatric Tracheostomy Tube (PDC) and the Cuffed Long Paediatric Tracheostomy Tube (PLC) feature gentle cylindrical cuffs.



#### ORDERING INFORMATION:

REF	Size	I.D. (mm)	O.D. (mm)	Length (mm)
<b>NEO</b>				
3.0NEO	3,0	3,0	4,5	30
3.5NEO	3,5	3,5	5,2	32
4.0NEO	4,0	4,0	5,9	34
4.5NEO	4,5	4,5	6,5	36
<b>PED</b>				
3.0PED	3,0	3,0	4,5	39
3.5PED	3,5	3,5	5,2	40
4.0PED	4,0	4,0	5,9	41
4.5PED	4,5	4,5	6,5	42
5.0PED	5,0	5,0	7,1	44
5.5PED	5,5	5,5	7,7	45
<b>PDL</b>				
5.0PDL	5,0	5,0	7,1	50
5.5PDL	5,5	5,5	7,7	52
6.0PDL	6,0	6,0	8,3	54
6.5PDL	6,5	6,5	9,0	56

Carton Quantity: 1

Shiley™ tracheostomy tubes are single packed and sterile.

#### ORDERING INFORMATION:

REF	Size	I.D. (mm)	O.D. (mm)	Length (mm)	Cuff Ø (mm)
<b>PDC</b>					
4.0PDC	4,0	4,0	5,9	41	12.4
4.5PDC	4,5	4,5	6,5	42	12.5
5.0PDC	5,0	5,0	7,1	44	15.0
5.5PDC	5,5	5,5	7,7	46	14.9
<b>PLC</b>					
5.0PLC	5,0	5,0	7,1	50	15.0
5.5PLC	5,5	5,5	7,7	52	14.9
6.0PLC	6,0	6,0	8,3	54	19.3
6.5PLC	6,5	6,5	9,0	56	19.3

Carton Quantity: 1