



**MAKE YOUR FIRST ATTEMPT
YOUR BEST.
EVEN FOR YOUR SMALLER PATIENTS.**

McGRATH™ MAC video laryngoscope for pediatrics

Medtronic
Further. Together

Pediatric airway complications are often unpredictable. And sometimes tracheal intubations (TI) require multiple attempts — especially if your patients are younger and have a history of difficult airway.¹ But you know that each additional attempt is associated with increased desaturation and TI-associated events (TIAEs).¹

You can be better prepared for unexpected airway complications with the McGRATH™ MAC video laryngoscope, which offers:

- Confidence — familiar Macintosh technique; portrait display video
- Simplicity — single-button functionality; minute-by-minute battery indication
- Convenience — disposable, sterile blades

In a study performed on pediatric mannequins, the McGRATH™ MAC video laryngoscope decreased intubation time and improved first time success rates over traditional Macintosh laryngoscopes.²

ORDERING INFORMATION

| SKU NO. | DESCRIPTION | SIZE | QUANTITY |
|-------------|-----------------------------|------|----------|
| 301-000-000 | McGRATH™ MAC handle | N/A | 1 |
| 350-072-000 | McGRATH™ MAC blade size 1 | 1 | 50/box |
| 350-017-000 | McGRATH™ MAC blade size 2 | 2 | 50/box |
| 350-005-000 | McGRATH™ MAC blade size 3 | 3 | 50/box |
| X3-003-000 | McGRATH™ MAC X blade size 3 | 3 | 10/box |
| 350-013-000 | McGRATH™ MAC blade size 4 | 4 | 50/box |
| 340-000-000 | 250-minute battery pack | N/A | 1 |

TECHNICAL SPECIFICATIONS

McGRATH™ MAC handle

| | |
|--------------|---|
| Size | 180 x 68 x 110 mm |
| Weight | 180 g |
| Power | Proprietary 3.6V lithium battery pack (c. 250 minutes) with auto shut off |
| Light source | High-intensity LED |
| Display | 2.5 in. LCD color display |
| Camera | CMOS |
| Materials | Durable medical-grade thermoplastics composite with structural alloy core. The device and packaging are latex free. |
| Warranty | 36 months |

McGRATH™ MAC disposable blades

| | |
|-----------|--|
| Material | Fog-free medical-grade optical polymer |
| Packaging | Packaged sterile for single use |

1. Lee JH, Turner DA, Kamat P, et al. The number of tracheal intubation attempts matters! A prospective multiinstitutional pediatric observational study. *BMC Pediatrics*. 2016;16:58.

2. Szarpak L, Truszczyński Z, Czyżewski L, Gaszyński T, Rodríguez-Núñez A. A comparison of the McGRATH-MAC and Macintosh laryngoscopes for child tracheal intubation during resuscitation by paramedics. A randomized, crossover, manikin study. *Am J Emerg Med*. 2016 Aug;34(8):1338-41.

© 2017, 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 11/2019-17-AW-0007(2)-[WF#2927634]



McGRATH™ MAC blade size 1



McGRATH™ MAC blade size 2



McGRATH™ MAC blade size 3